

## CXCL15/ Lungkine (C-6His), Mouse, Recombinant

货号 : PCK235

### 产品信息

别名	C-X-C Motif Chemokine 15; Cxcl15; Lungkine; Small-inducible Cytokine B15; Scyb15
物种	Mouse
表达宿主	Human Cells
序列信息	Gln26-Ala167
检索号	Q9WVL7
分子量	17.4 kDa
标签	C-6His

### 产品特性

纯度	>95% as determined by reducing SDS-PAGE.
内毒素	< 1.0 EU per µg as determined by LAL test.
保存	Lyophilized protein should be stored at -5~-20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -5~-20°C for 3 months.
运输	Ambient temperature or ice pack.
制剂	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.



## 复融

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## 背景介绍

Mouse C-X-C Motif Chemokine 15, also known as Cxcl15, is a secreted Protein which is member of the ELR Motif-containing CXC Chemokines. It expressed at low levels in fetal lung, the expression restricted to the lung, produced by bronchoepithelial cells and is released into the airways. It plays an important role in lung-specific neutrophil trafficking during normal and inflammatory conditions. It also appears chemotactic for neutrophils.

## SDS-PAGE

